

The Bootie Butler® Shoe Cover System makes dispensing shoe covers a simple, safer and money-saving process.

Dispensing shoe covers is now as simple as 1-2-3:

Dispensar las cubiertas de zapatos ahora es tan simple como 1-2-3:

01/



Place one foot in the dispenser, heel first.

Coloque un pie en el dispensador, primero el talón.

02/



Slowly pull back until cover stretches beyond shoe tip, point toe down and pull back.

Lentamente tire hacia atrás hasta que la cubierta se extienda más allá de la punta del zapato, apunte hacia abajo y tire hacia atrás.

03/



Repeat with other foot.

Repita con el otro pie.